

Prosthodontic

Wax & Modelling Carvers

Cement Spatulas

Wax Knives & Spatulas

Impression Trays

Crown Removers




WAX AND MODELLING CARVERS




Zahle WCZ1	Zahle WCZ2	Zahle WCZ3	Hylin WCH1	Hylin WCH2	Hylin WCH3	Gritman WCG	Lecron WCL1	Lecron WCL4
---------------	---------------	---------------	---------------	---------------	---------------	----------------	----------------	----------------

PKT CARVERS COLOUR CODED



PKT CARVERS



WAX KNIVES & SPATULAS



Fahnenstock
WKFS



Lessmann
WKLS



Fahnenstock
WKFL



Lessmann
WKLL



Knife
OK5A



Gritman
OKG



Flexible
PS6




Spatula
PS3



Spatula
PS4

HEIDEMANN SPATULA



PFIHS1 #1 #6
 PFIHS2 #1 #6
 PFIHS3 #1 #6

CEMENT SPATULA



CS6 #3
 CS5 #3
 CS4 #3 #4
 CS3 #3

GINGIVAL CORD PACKER



Serrated GCP113 #3 #4
 Non-Serrated GCP113N #3 #4
 Serrated GCPAS113 #3 #4
 Serrated GCPR7 #3 #4
 Non-Serrated GCPR7N #3 #4

IMPRESSION TRAYS DENTULOUS PERFORATED



DENTULOUS NON-PERFORATED



IMPRESSION TRAYS EDENTULOUS PERFORATED

UPPER



ITRLEDPU1

ITRLEDPU2

ITRLEDPU3

ITRLEDPU4

LOWER



ITRLEDPL1

ITRLEDPL2

ITRLEDPL3

ITRLEDPL4

EDENTULOUS NON-PERFORATED

UPPER



ITRLEDNPU1

ITRLEDNPU2

ITRLEDNPU3

ITRLEDNPU4

LOWER




ITRLEDNPL1

ITRLEDNPL2

ITRLEDNPL3

ITRLEDNPL4

IMPRESSION TRAY DENTULOUS PERFORATED PEEDO



IMPRESSION TRAYS HALF



ANTERIOR AND POSTERIOR



ITRLDPUA



ITRLDPLA

ADJUSTABLE



ITRLDPA



ITRLDNPA

CALIPER IWANSON



Metal
CLP1



Wax
CLP2

ARTICULATING PAPER HOLDER



Miller Straight
APF2E



Miller Curved
APF2C

CROWN REMOVER

TURKEY PATTERN



CRT

CRT1

CRT2

CRT3

CRT4

CRT5

CRT6

CRT8

CRT9

MANUAL PATTERN



CRMP

CRMP1

CRMP2

CRMP3

CROWN REMOVER

AUTOMATIC



CRS

CRP

Crown Remover Tips with CRS & CRP



CR1

CR2

CR3

CR4

CR5

CR6

CR7
(For CRS)

CR8
(For CRP)

CR9

CR10

Basic Kit Includes :

- Crown and Bridge Remover # CRS/CRP
- Single crown remover tip # CR1
- Crown remover tip Bifid # CR2
- Crown & bridge removing tip Sickle # CR3
- Single crown remover Contra angle tip # CR4
- Single crown remover tip Inverted curved # CR5
- Single Crown remover tip straight # CR6
- Bridge remover Flexi-wire set # CR7 (For CRS)
- Bridge remover tip with lock # CR8 (For CRP)
- Tips adjustment Key # CR9
- Open-end adjustment wrench # CR10

Terms of Use

Removal of :

- Permanently cemented crowns and /or bridges
- Temporarily cemented provisional or final restorations
- Maryland bridges
- Orthodontic brackets/bands
- Removable appliances "stuck" at attachments
- Multi-unit restoration "stuck" on working lab model

Precautions :

1. To remove crowns placed over pins or screw type posts, use lowest "strike intensity" to prevent fracture.
2. This instrument could cause damage to the margins of soft gold crowns. If re-use is anticipated, cut a small groove parallel to the lingual margin.
3. Highest strike intensity could cause root fracture; be careful in long span restorations.
4. Only low intensity strike force should be used with teeth having high mobility and little bone support.
5. For highest strike intensity settings, let the individual situation and common sense prevail.

ARTICULATOR



A5222

CROWN HOLDING FORCEP



CHFRS
12.5cm

CHFRC
12.5cm

PUTTY CUTTER



Right
Handress
PCR

Left
Handress
PCL

PLASTER SHEARS



PCS
21cm

CROWN OPENER



CR4566

TEMPORARY
CROWN REMOVING PLIER



CR4562